


 **REPRODUCTIVE GENETIC INNOVATIONS LLC**


 **PGDIS**

IMPACT of 24-chromosome ANEUPLOIDY TESTING on the OUTCOME of PGD for MONOGENIC DISORDERS

RECHITSKY S, KULIEV A.



 **REPRODUCTIVE GENETIC INNOVATIONS LLC**

 **PGDIS**


APPROACHES TO Universal PGD

Combined PGD :

- **SGD and Aneuploidy on the same sample**

NOW : Combination of Different types of testing

FUTURE : One method for all conditions





PRESENT RGI EXPERIENCE



TOTAL : Over 12,500 CASES

4501 PGD CYCLES

FOR 426 MENDELIAN CONDITIONS

1479 Single Gene cases + Aneuploidy

762 Single Gene cases + 24 AT

< 900 TRANSLOCATIONS



FLOW CHART FOR COMBINED TESTING



A single cell is placed in lysis buffer and spanned

Mutation & Linkage (Haplotyping):
PCR Mix with multiple outside primers,

Second Round started as a separate PCR reaction for each locus





Combined testing: SGD+AT-24

WGA:
aCGH+Mutation --: PEP,DOP





Karyomapping : MDA ----->

NGS: DOP,MDA or MALBAC ----->



Single Gene Analysis:

Aneuploidy testing





WGA APPROACHES :

- 1. PCR BASED: **PEP-PCR** (15 MER DEGENERATE PRIMER)
DOP-PCR (QUASI-RANDOM ANNEALING OF HYBRID OLIGOS)
LA-PCR (ADAPTORS WITH UNIVERSAL SEQUENCE LIGATED TO DNA ENDS)
- 2. ISOTHERMAL: **MDA** (RANDOM PRIMING AND STRAND DISPLACEMENT)
- 1. COMBINATION OF BOTH: **MALBAC** (DISPLACEMENT PRE-AMPLIFICATION and LOOP CREATION)


MICROARRAY or NGS limitations due to WGA nature:



- BREATH OF GENOMIC COVERAGE;
- AMPLIFICATION BIAS DUE TO %GC CONTENT;
- ALLELE DROP OUT (ADO)
- PREFERENTIAL ALLIIC AMPLIFICATION (PA)
- CHIMERIC DNA MOLECULES;
- NUCLEOTIDE COPY ERRORS
- CAN VARY SIGNIFICANTLY FOR DIFFERENT WGA APPROACHES AND CAN AFFECT INTERPRETATION


A ONE-FIT- ALL WGA METHOD REMAINS TO BE ESTABLISHED



ADO RATE IN DIFFERENT TYPE OF CELLS with or without WGA OR MDA















Test type	Total # families	Relatives were Not available for haplotype prediction	Same mutation (including offspring)	De novo mutations	Recombinant baby
SGD	2460	354 12.73%	41 1.7%	179 7.2%	8 0.3%
HLA Typing	320	NA	NA	NA	17 5.3%
TOTAL	2780	354 (12.73)	41 (1.7%)	179 (7.2%)	25 (0.89%)



Total : 22.52%











Overall Results and outcome of PGD for Single Gene Disorders, HLA & Aneuploidy Testing

Testing	Patient/ Cycle	# ET	# Embryos	Pregnancy	SAB	Delivery	# Babies
SGD (INCLUDING HLA)	1749/2791	2113	3881 1.8	1078 51.0%	158 14.6%	920	1150
SGD+A by PCR or FISH	492/948	707	1209 1.7	323 45.6%	30 9.2%	293	357
SGD+ 24 AT For 122 conditions (INCLUDING HLA)	485/762	463*	616 1.3	376 81.2% <small>P<0.0001</small>	29 7.7% <small>P<0.0027</small>	347 92.3%	382
TOTAL	2726/4501	3283	5706	1777	217	1560	1889

CONCLUSIONS:

- **Trophectoderm Biopsy** (Day 5&6) has improved amplification efficiency (fewer inconclusive samples) and provide more embryos suitable for transfer
- **WGA** allowed multiple conditions testing on the same sample and embryo identification
- **Adding aneuploidy testing by aCGH or NGS**
 - Increased pregnancy rates
 - Reduced chance of miscarriage
 - Predicted accuracy is 98%



RGI TEAM

PREIMPLANTATION MOLECULAR GENETICS

Svetlana Rechitsky, Ph.D.
Oleg Verlinsky, M.A.
Tatiana Pakhalchuk B.S
Adam Goodman, B.S
Geraldine San Ramon ,M.S
Maria Prokhorovich, PhD
Rita Gershman, B.S
Michael Ducor, M.S
Lukas Machaj, B.S.

CLINICAL GENETICS

JL Simpson, M.D., FACMG
Lee Shulman, MD., FACMG
Dyvia Shah, M.S.
Agness Machaj, B.S., M.S.
John Zimmerman, B.S.

PREIMPLANTATION CYTOGENETICS

Anver Kuliev ,MD, PhD
Yuri Ilkevitch, PhD.
Zev Zlatopolsky M.S.

ASSISTED REPRODUCTION TECHNOLOGY AND MICROMANIPULATION

Svetlana Lerner M.S., Natalie Ilkevitch M.S.
Yury Ilkevich, PhD; Nathan Katz, PhD